

# BI@ BS3

## Portable Diagnostic Dermatoscope

### Product specifications:

Model: BS3

Product ratio: 10 times

Product Scale: 2cm accuracy: 0.5mm / grid

Weight: 225 g

Product Size: 180mm (cylinder head Length), 41mm (diameter), (55mm bottom diameter)

Product Material: Metal, Optical glass lens

Battery Type: 2\* AA batteries (without battery factory)

Packing: cloth cover, white box

Packing materials: 10 units / inner box, 40 units / box,

Carton weight: 42.5 \* 34.5 \* 36.5CM, 12 / 9KG



### Product Operation:

- 1 install AA batteries in the BS3 digital microscope Lane;
- 2 hold the power button for two seconds to on;
- 3.observe the objects will be displayed on the screen.
- 4.adjust the ratio until the draw the clearest image of the object; in clearly, And let the image alienment the calibration.
- 5.turn off the power button and put in back to the box.



**Wuxi Biomedical Technology Co., Ltd.**

[www.portable-ultrasound-scanner.com](http://www.portable-ultrasound-scanner.com)

Add:Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel:86-510-85572789

Fax:86-510-85572789

Email:[sales@cnbiomedical.com](mailto:sales@cnbiomedical.com)

Mobilephone:86-13382885910



is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.